



Multi-parameter Indicator Class 4



SAP / Monitoring



Multi-Parameter Indicator

REF: 104.003.0250 Multi-parameter indicator 134°C-3,5 min Unit: 250 pieces

The multi-parameter indicator is checking time and temperature of the sterilization cycle at the exact place where the indicator is placed. A clear blue to black colour change is indicating the correct steam sterilization cycle parameters have been met.

The Sterintech™ Indicators are intended to be used during steam sterilization to check critical parameters at the place where they are placed.

Characteristics:

- Shelf life : 5 years after production
- Packed in an easy dispensing carton
- Steam indicator conform to ISO 11140-1 : Class 4
- Clear color transition into black
- Chemical indicators covered with lacquer
- Non-Toxic / Lead Free
- Other types and lay-outs available upon request.

