



# Self-Adhesive Indicator Strip



## SAP / Monitoring

<p>STEAM Indicator ISO 11140-1 - class 4</p> <p>Green = Sterilized REF : 104.010.0500</p>	<p>STEAM Indicator ISO 11140-1 - class 4</p> <p>Green = Sterilized REF : 104.010.0500</p>
<p>Formaldehyde Indicator ISO 11140-1 - class 4</p> <p>Green = Sterilized REF : 104.040.0500</p>	<p>Formaldehyde Indicator ISO 11140-1 - class 4</p> <p>Green = Sterilized REF : 104.040.0500</p>
<p>Ethylene Oxide Indicator ISO 11140-1 - class 4</p> <p>Brown = Sterilized REF : 104.030.0500</p>	<p>Ethylene Oxide Indicator ISO 11140-1 - class 4</p> <p>Brown = Sterilized REF : 104.030.0500</p>
<p>Plasma Sterilization Indicator acc. ISO 11140-1</p> <p>Yellow = Sterilized REF : 104.020.0500</p>	<p>Plasma Sterilization Indicator acc. ISO 11140-1</p> <p>Yellow = Sterilized REF : 104.020.0500</p>

## If you want to document .....

REF: 104.010.0500 STEAM Sterilization indicator	- Class 4 Unit: 500 pieces / roll
REF: 104.021.0400 H <sub>2</sub> O <sub>2</sub> Sterilization indicator	- Unit: 400 pieces / roll
REF: 104.030.0500 EtO Sterilization indicator	- Class 4 Unit: 500 pieces / roll
REF: 104.040.0500 FORMALDEHYDE Sterilization indicator	- Class 4 Unit: 500 pieces / roll

These indicator strips can be used as in-pack indicator in order to detect sterilization conditions have been reached at any local position. They come as a piggy-back self-adhesive label in order to document these into the patient records.

### Characteristics:

- Shelf life : 3 years after production
- On easy dispensing roll , Dimensions: 70 x 20 mm (per strip)
- Conforms to ISO 11140 - 1 (Class 4 for all except the Plasma version )
- Very easy interpretation of result
- Complete indicator area sealed by lacquer
- Non-Toxic / Lead Free

