



Sequential Emulator Class 6



SAP / Monitoring



Emulator with a big plus

REF: 106.003.0250 Sequential Emulator 134°C - 3,5 min Unit: 250 pieces
REF: 106.007.0250 Sequential Emulator 134°C - 7,0 min Unit: 250 pieces

Emulators can be used in addition to the batch control Helix which will allow you to release load. Especially when OR staff need to be re-affirmed about the sterilization conditions.

The Sterintech™ Emulators are offering you more information about the steam sterilization process than those with only one emulator printed. Due to that the three emulators are calibrated to different time/temperature combination it is providing information about for instance over-heating.

Characteristics:

- Shelf life: 5 years after production
- Packed in an easy dispensing carton
- Steam emulator conform to ISO 11140-1 : Class 6
- Clear colour transition into black of the emulator
- Chemical emulator covered with lacquer
- Non-Toxic / Lead Free

