



Traceability Two Line Label Gun



Documentation



Two line but it carries all the information

Two lines ... but it carries all information

The Sterintech two line label gun is reliable, easy to use and durable. It will serve for years of intensive usage. Sterintech double labels for two line label guns are available with all type of sterilization indicators as well as a variety of colors. Colored labels are used to assist FIFO systems as well as easy recognition of sterile goods for a certain department.

Often the combination of a double label gun label and normal tape is a more cost-effective solution than steam tapes only. Next to that they offer better readable data and prevent writing errors.

Characteristics:

- Process Indicators (Steam, FO, EtO, Plasma) conform to ISO 11140 – 1 Class 1
- Clear transition of indicator into its end-color
- Easy to use label gun
- Double labels serve for an easy traceability and documentation
- Non-Toxic / Lead Free indicators
- Latex free label materials
- Can be fully recycled

For reference numbers and types please refer to the next page.





Traceability Two Line Label Gun



Documentation



Simplicity of data recording

Steritech two line label gun:

- REF: 602.003.0001 Two line label gun Unit: 1 piece
- REF: 609 102 0001 Two line label gun Ink roller Unit: 1 piece

Steritech Two Line labels are standard 36mm x 19mm and available as follows:

- REF: 602.011.9000 Label with Steam indicator Violet-Green Unit: 9.000 labels
- REF: 602.021.6000 Label with Plasma indicator Red-Yellow Unit: 6.000 labels
- REF: 602.030.9000 Label with EtO indicator Yellow-Brown Unit: 9.000 labels
- REF: 602.040.9000 Label with Formaldehyde indicator Violet-Green Unit: 9.000 labels
- REF: 602.101.9000 Label with Blue colour Unit: 9.000 labels
- REF: 602.102.9000 Label with Red colour Unit: 9.000 labels
- REF: 602.103.9000 Label with Green colour Unit: 9.000 labels
- REF: 602.104.9000 Label with Yellow colour Unit: 9.000 labels
- REF: 602.105.9000 Label with Orange colour Unit: 9.000 labels
- REF: 602.106.9000 Label with Grey colour Unit: 9.000 labels
- REF: 602.107.9000 Label with Purple colour Unit: 9.000 labels
- REF: 602.114.9000 Label with Steam indicator Violet-Green & Yellow colour Unit: 9.000 labels

