



Container Filter



Packaging



Easy to use Container Filters

The Sterintech Container Filter is based and compliant to ISO 11607 and ISO 11140 standards. The used materials are mechanically strong to withstand the vacuum and steam pulse movements and is an effective biological barrier. Available in two versions:

- REF: 901.101.0500 Container Filter round central indicator Unit: 500 pieces
- REF: 901.102.0500 Container Filter round side indicator Unit: 500 pieces

Characteristics:

- Medical grade paper , 60 gram per square meter.
- In accordance with ISO 11607 - 1 : Hydrophobic, for containers
- High mechanical strength and low liquid permeability
- Steam indicator conform to ISO 11140-1 : Class 1
- Clear colour transition of the indicator (Yellow - Dark)
- Contains no toxic materials
- Other types available upon request.



Steam Spots:

- REF: 904 001 0010 Steam Spots 16 x 9mm Din A4 Sheets with 340 spots Unit: 10 pieces