



Indicator Tape



Packaging



Superior Quality / Super Stick

To maintain a constant quality at all production stages all processes are subject to ISO 13485 procedures. No need to mention that our steam tapes are tested and certified per batch before distributing them to our customers.

| | | |
|---|----------------|-----------------|
| REF: 502.019.0048 Steam Indicator Tape | 19 mm x 50 mtr | Unit: 48 pieces |
| REF: 502.025.0036 Steam Indicator Tape | 25 mm x 50 mtr | Unit: 36 pieces |
| REF: 503.019.0048 EtO Indicator Tape 600mg/Ltr | 19 mm x 50 mtr | Unit: 48 pieces |
| REF: 504.025.0036 H ₂ O ₂ Plasma Indicator Tape | 25 mm x 50 mtr | Unit: 1 piece |

Providing a quick recognition of sterile goods by its indicator while securing the sterility of wrapped items our steam tape performs in Hospitals anywhere in the world.

Characteristics:

- Single packed to guarantee shelf-life specifications up to 2 years after production
- Easy unwinding from the roll
- Steam indicator conform to ISO 11140 - Class 1
- Steam indicator conform to ISO 13485
- Latex free glue
- Non-Toxic / Lead Free

