



## Sterilization Rolls



### Packaging



### Sterile barrier confidence .....

Today Infection control is getting increasingly important since more and more bacteria starts to get resistant to the latest anti-biotic's. Infection Prevention is the key word. Sterilization packaging materials have been developed over many years to increase the stability of the sterile barrier over time.

Sterintech™ Sterilization rolls are used to wrap medical devices. The flexibility of length cutting is making a dedicated packaging of individual instruments possible and offer a number of advantages to the end-user

Both flat and gusseted sterilization rolls are offering two indicators for STEAM and EtO. Quality and selection of the paper and film is providing a ISO 11607 compliant sterile barrier.

The Sterintech™ Sterilization rolls are characterized by the following:

- Easy recognition of content by the multi-layer transparent film
- Strong, multi-line seals
- Indicators and text printed outside the packing area
- Non-tearing film avoids tearing at opening the pack
- Clean and fire-free opening of pack
- Compatible to STEAM, FORM or EtO sterilization by three imprinted indicators
- Compliant to standard ISO 11607-1, EN 868 - 2 and ISO 11140
- Please refer to the next page for all standard dimensions.





# Sterilization Rolls



## Packaging



### Sterile barrier confidence .....

Standard dimensions Flat reel:

- REF: 301.050.0001 Flat reel, 60 gr/sqm, 50 mm x 200 mtr Unit: 1 pcs
- REF: 301.075.0001 Flat reel, 60 gr/sqm, 75 mm x 200 mtr Unit: 1 pcs
- REF: 301.100.0001 Flat reel, 60 gr/sqm, 100 mm x 200 mtr Unit: 1 pcs
- REF: 301.150.0001 Flat reel, 60 gr/sqm, 150 mm x 200 mtr Unit: 1 pcs
- REF: 301.200.0001 Flat reel, 60 gr/sqm, 200 mm x 200 mtr Unit: 1 pcs
- REF: 301.250.0001 Flat reel, 60 gr/sqm, 250 mm x 200 mtr Unit: 1 pcs
- REF: 301.300.0001 Flat reel, 60 gr/sqm, 300 mm x 200 mtr Unit: 1 pcs
- REF: 301.350.0001 Flat reel, 60 gr/sqm, 350 mm x 200 mtr Unit: 1 pcs
- REF: 301.400.0001 Flat reel, 60 gr/sqm, 400 mm x 200 mtr Unit: 1 pcs

Standard dimensions Gusseted reel:

- REF: 302.075.0001 Gusseted reel, 60 gr/sqm, 75 mm x 100 mtr Unit: 1 pcs
- REF: 302.100.0001 Gusseted reel, 60 gr/sqm, 100 mm x 100 mtr Unit: 1 pcs
- REF: 302.150.0001 Gusseted reel, 60 gr/sqm, 150 mm x 100 mtr Unit: 1 pcs
- REF: 302.200.0001 Gusseted reel, 60 gr/sqm, 200 mm x 100 mtr Unit: 1 pcs
- REF: 302.250.0001 Gusseted reel, 60 gr/sqm, 250 mm x 100 mtr Unit: 1 pcs
- REF: 302.300.0001 Gusseted reel, 60 gr/sqm, 300 mm x 100 mtr Unit: 1 pcs
- REF: 302.350.0001 Gusseted reel, 60 gr/sqm, 350 mm x 100 mtr Unit: 1 pcs
- REF: 302.400.0001 Gusseted reel, 60 gr/sqm, 400 mm x 100 mtr Unit: 1 pcs

Bespoke products available upon request.

