



Tyvek Sterilization Rolls



Packaging



Standard size Tyvek Rolls for H₂O₂ Plasma Sterilizers

Tyvek has been introduced into the Medical Device Industry as packaging and nowadays finds its place at CSSD departments where H₂O₂ Plasma sterilization is taking place.

The Sterintech™ Tyvek Sterilization rolls are based upon a 74 gr/sqm Tyvek material in combination with a 52 µm strong PE / PET film providing a secure sterile barrier during and after Plasma sterilization.

Available in a big variety of widths varying from 75 to 400 mm the Tyvek Sterilization rolls come in 70 mtr rolls.

The Sterintech™ Tyvek Sterilization rolls are characterized by the following:

- Easy recognition of content by the multi-layer transparent film
- Indicator and text printed outside the packing area
- Non-tearing film avoids tearing at opening the pack
- Suitable for H₂O₂ Plasma and EtO sterilization (only plasma indicator provided)
- Compliant to standard ISO 11607-1, EN 868 - 2 and ISO 11140

(Tyvek is a trademark of Dupont)





Tyvek Sterilization Rolls



Packaging



Standard size Tyvek Rolls for H₂O₂ Plasma Sterilizers

Standard Tyvek reel dimensions for H₂O₂ Plasma Sterilizers:

- REF: 303.075.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 75 mm x 70 mtr Unit: 1 pcs
- REF: 303.100.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 100 mm x 70 mtr Unit: 1 pcs
- REF: 303.150.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 150 mm x 70 mtr Unit: 1 pcs
- REF: 303.200.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 200 mm x 70 mtr Unit: 1 pcs
- REF: 303.250.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 250 mm x 70 mtr Unit: 1 pcs
- REF: 303.300.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 300 mm x 70 mtr Unit: 1 pcs
- REF: 303.350.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 350 mm x 70 mtr Unit: 1 pcs
- REF: 303.400.0001 Tyvek reel, 74 gr/sqm, 52µm PET/PE, 400 mm x 70 mtr Unit: 1 pcs

(Tyvek is a trademark of Dupont)

