



## Bowie Dick Helix Test (SS)



### SAP / Monitoring



## Checks your sterilizer for various loads .....

REF: 101.003.0100 SS Helix with 4 mtr tube Unit: 100 indicators  
REF: 101.004.0250 SS Helix with 4 mtr tube Unit: 250 indicators

The Sterintech Bowie Dick Helix Test with Stainless Steel Helix and 4 mtr long tube can be used as Bowie Dick test for both load types: Hollow and Porous Loads.

Helix devices and their indicators have proven that they detect failures during each cycle. From leakages, bad steam conditions, non-condensable gases, time and temperature issues and forms a simple aid to release your load securely.

### Characteristics:

- Conform to ISO 11140-1 Class 2 and ISO 11140-4 : Class 2
- Clear transition of colour into black
- Easy to interpret
- Easy to document
- Non-Toxic / Lead Free

