



Bowie Dick Test Sheet



SAP / Monitoring



Once a day monitoring own pack

REF: 102.010.0050 Unit: 50 pieces

The Steritech Bowie Dick Test Sheet with its DIN A4 dimension is designed for end-users who are creating their own Bowie Dick Test Pack. The laquer laminated indicator sheet is able to be filled in with all relevant data.

Used for daily monitoring of pre-vacuum steam sterilizer. Based upon the Bowie Dick test result the steam sterilizer is released for usage during the day.

Characteristics:

- To be packed into pack as per standard / instructions
- Fulfil requirements of EN 285
- Conform to ISO 11140-1 and ISO 11140-3: Class 2'
- Dimension: DIN A4 size
- Clear yellow – black colour change
- Easy to interpret
- Easy to document
- Non-Toxic / Lead Free

The Bowie Dick Test Sheet detects failures like:

- Air leaks,
- Insufficient Vacuum,
- Inadequate Steam Penetration.