



EtO / Formaldehyde Helix Test



SAP / Monitoring



Helix Tests for EtO and Formaldehyde Sterilizers

REF: 103.030.0100 Hollow A Helix with 1,5 mtr tube Unit: 100 indicators (EtO)

REF: 103.040.0100 Hollow A Helix with 1,5 mtr tube Unit: 100 indicators (Formaldehyde)

Process Challenging Device with a 1,5 mtr Teflon tube, which is used to test the cycle performance of EtO and Formaldehyde.

The Helix Test for EtO and Formaldehyde is able to detect:

- Sterilization effectiveness too low
- Sterilization effectiveness too short
- Insufficient vacuum in depth and in number of vacuum pulses.
- Insufficient air removal from hollow A devices
- Insufficient penetration in hollow devices
- Leakage of piping / valves / door seals
- Detection of presence of small volume inert gases (or air)

Characteristics:

- Conform to ISO 11140-1 Class 5
- Clear transition of colour
- Easy to interpret
- Easy to document
- Non-Toxic / Lead Free
- Can be fully recycled

