



MEDI-CLAVE Pty (Ltd)
Turnkey solutions for today's cssd

Head Office
Portion 69, Melvic
Schurberg Ext 488
Pretoria, South Africa

PO Box 1652
Hartbeespoort
0216, South Africa
Tel:+27(0)12 371 3759
Fax:+27 (0)12 551 7928
admin@medi-clave.co.za
sterintech@medi-clave.co.za
Website: www.medi-clave.co.za
Registration No. 1999/018374/07

Protein Rapid Swab Test

Monitoring



Curious about your protein level ...

The EN-ISO-15883 standard is recommending a routine performance check of your washer disinfector. Rather than checking the cleanliness of your instruments visually you can easily test it by checking the level of proteins. The Protein Rapid Swab Test allows you to check the level of proteins left after washing disinfection on your instruments without any laboratory costs and within 5 seconds. The yellow swab will turn green when proteins are detected above 20 μ g

REF: 906.003.0050!!Protein Rapid Swab Test! !Unit: 50 double swabs (100 tests)

Characteristics:

- ◆ Detects 20 μ g of proteins, well below the required 100 μ g as per EN ISO 15883
- ◆ Easy to use, easy to interpret, no training required
- ◆ Very cost-effective in relation to laboratory tests
- ◆ No guesswork anymore by visual check but knowing for sure
- ◆ Quick test results: less than 5 seconds
- ◆ Non-toxic - safe to use - easy disposal.
- ◆ Calibrated test results
- ◆ Shelf-Life: 2 years after production date



 **Sterintech**™

Sterilizing is our Business! and our only Business!

