



Container Seal



Packaging



Anti-Tamper Container Seal with Process Indicator

Container Seals:

In order to secure the container for opening by un-authorized staff the container seals are offering an important shackle in the quality management chain.

Prevents opening by un-authorized staff. Only by breaking the seal the container can be opened. A check at the OR for broken seals is guaranteeing the sterility of materials from the moment of sterilization till the actual use. The Sterintech™ Container seals are the only ones available with a tested ISO 11140 class 1 Steam process indicator.

REF: 902.101.1000	Container seal Plastic Blue	Unit: 1000 pieces
REF: 902.102.1000	Container seal Plastic Blue with Indicator	Unit: 1000 pieces
REF: 902.106.1100 SA	Container seal Paper with Indicator & Pin 29x65mm&140mm pin	Unit: 1100 pieces

Characteristics:

- Produced in polypropylene , softening above 150°C.
- Suitable for steam sterilization processes.
- Fulfilling EN 868 part 8 section 4.2.3
- Can be fully recycled.
- Steam indicator conform to ISO 11140-1 : Class 1
- Clear color transition of the indicator
- Other types and colors available upon request.

